1C / 5 / UDF
CON UGELLO A DOPPIO FLUIDO
1C / 5 / UDF-T
CON UGELLO A DOPPIO FLUIDO E TURBINA

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I.C.F. LABORATORY SPRAY DRYER has been conceived to supply the technological and research laboratories with an instrument having a duty being as much flexible as possible, in order to examine the behaviour of products and formulations in solution, emulsion or suspension, during the spray drying process.

OPERATION PRINCIPLE

The liquid product is "spray dried", i.e. reduced to drops being as small so as to obtain a very wide specific surface. The drying air acts on so wide a surface that it allows an instantaneous evaporation of the liquid part with a consequent immediate fall in air temperature. Thanks to this principle, it is possible to process in the spray dryer thermolabile products without any loss of activity.

ATOMIZING SYSTEMS AND DESCRIPTION OF THE PROCESS.

The atomizing systems of the laboratory spray dryer are two and well distinct:

- disk or centrifugal rotor in co-current configuration

- double-fluid nozzle (product and air or other gases) in co- or counter - current configuration.

It is extremely easy to change the number of disk revolutions, and the pressure of atomizing compressed air of the double-fluid nozzle. The liquid to be spray dried is batched by means of an adjustable-delivery pump. The atomized liquid comes into contact with a stream of hot air generated by gas or electric resistances. The dried product is collected inside a container placed below the drying chamber: the finest or surplus product is sucked together with the process air and collected below the high-efficiency cyclone separator. A centrifugal fan sucks the exhaust air and discharge it in the atmosphere. The laboratory spray dryer is also provided with an electric control board including all the apparatus for checking, protection, control and regulation.

TECHNICAL DATA

- Evaporating capacity: 3-6 kg/h of water
- Inlet-air temperature: • up to 330° with gas heater
• 220° with electric heater
- Outlet-air temperature: 60-150° C
- Installed electric power: • 1,3 kW with double fluid nozzle
• 2,8 kW with double fluid nozzle and turbine
- Heating thermic capacity: 10.000 Kcal/h with gas - 12 kW with electric heater
- Compressed-air consumption: 500 NI/min. at 6 bar max
- Chamber inside diameter: 1.000 mm.
- Max. width: 2.000 mm.
- Max. length: 2.000 mm.
- Max. height: 3.000 mm.

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I.C.F. & Welko S.p.A.